

Authorized

**PAN India Distributor & Importers**



**SANARE  
MEDICAL**



**PROGENY**  
Rehab tech

## (code CYB.SHOCK) Shock waves therapy

**PRO-SHOCK WAVES** it is a very advanced product delivering radial shock waves for therapeutical use and it is absolutely one of the most versatile device of the market.

The use of modern mechanical technologies applied to this machine guarantee the best quality of delivery of the shock wave on the patient. Furthermore, there are several energy transmitters with different focalizations, adapting the use of the machine to different types of application on the human body.

From sports injuries to chronic diseases and disorders requiring rehabilitation, the **PRO-SHOCK WAVES** offers the ultimate technology in shock waves and, with its wide range of pre-set protocols, provides highly effective results in the treatment of pain.

The unique patented **E.S.A.® System (Extremely Soft Absorbed)** reduce the energy transmitted to the therapist up to the 90%. With its innovative and patented System, the PRO-SHOCK WAVES provides unequalled parameter flexibility, allowing the therapist to reduce the energy and create the perfect combinations of parameters for each patient.



### SPECIFICATIONS

- Selection by anatomical area of the body
- Selection by medical branches
- Pulsed mode up to 25 Hz
- Pressure power up to 5 bar
- Possibility to create custom programs
- Patient Database
- Huge quantity of preset programs
- E.S.A.® System - This exclusive system reduce more than 90% of energy transmitted to the therapist during the treatment. This system works as a classic shock absorber and removes any problems from the transfer of energy from the handle to the therapist
- Diagram program



Equipment compliant with Directive 93/42/EEC integrated by Directive 2007/47/EC



### STANDARD EQUIPMENT

- 1 pc. - 15 mm focused transmitter [1] (A01.00142)
- 1 pc. - 15 mm multi-focused transmitter [2] (A01.00140)
- 1 pc. - 9 mm multi-focused transmitter [3] (A01.00141)
- 1 pc. - Gel tube [4] (R/273)
- 1 pc. - Gun case kit [5] (VAL-24/4)
- 1 pc. - Barrel Interchangeable Kit (A01.00104)
- 1 pc. - Key for locking ring nut [6] (J02.01717)



### OPTIONAL EQUIPMENT

- 36 mm transmitter for aesthetic medicine [7] (A01.00107)
- Trolley Model Z [8] (A01.00128)

### TECHNICAL DATA

Number of output:	1
Power supply:	230 Vac 50 Hz, ±10%
Maximum absorption:	150 W
Fuse protection on the mains:	230 Vac 2x2 A-T - 5x20 mm
Indications:	Backlight graphic LCD display 640 x 480 pixels Touch Screen
Type of emission of the pulses:	Single, Continuous, Burst
Type of emission Burst:	(1÷100) burst pulses (10÷2000) ms pause between burst pulses
Pressure:	1÷5 bar to step of 0.1 bar
Emission frequency:	1÷25 Hz to step of 1 Hz
N. of pulses emitted:	0÷10000
Weight:	8.85 Kg
Dimensions (LxWxH):	33 x 33 x 24 cm
Insulation class:	Class I - BF type



**SANARE MEDICAL PVT LTD**

Website : [www.sanaremedical.com](http://www.sanaremedical.com)

Email : [info@sanaremedical.com](mailto:info@sanaremedical.com)

Tel : +91 9810275300 , +91 98101 37486

Corp. Office: E-354 Greater kailash Part 1  
New Delhi 110048

Regd. Office : DP 61 Phase 1 Udyog Vihar  
Gurugram Haryana 122001