



HOSPITAL FURNITURE

www.sanaremedical.com



SANARE MEDICAL PVT LTD

**Features:**

- Made up of CRCA tube and CRCA Sheet.
- Available with four section 40 density form mattress
- Detachable/ Mounting rod polymer moulded head and foot panel
- Positions: Back rest, knee rest, trendelenburg/reverse trendelenburg, high – low, cardiac chair position and one touch zero level (CPR) operated by engraved nurse control panels on both panel
- Provided with 4 sectional polymer moulded safety side rails with separate angle indicators.
- Corner PU buffers absorb impact and reduce damage in handling.
- Provided with heavy duty saline stand with provision for fixing it at all four locations.
- Castors of 125/150 mm dia. with diagonally locking mechanism
- Dual side manual CPR
- Frame size (approx): 2030L X 1060 W X 450-750 H
- Back rest tilt : 0-75 degree
- Knee rest tilt : 0-35 degree
- Trendelenburg : 0-12 degree
- Reverse trendelenburg : 0-12 degree
- The complete metal components are pre treated and epoxy powder coated.

Optional: Central Locking System



Feature

s:

- Made up of CRCA tube and CRCA Sheet.
- Available with four section 40 density form mattress
- Castors of 125/150 mm dia. With diagonally locking mechanism.
- Position: Back rest, knee rest, trendelenburg / reverse trendelenburg, high – low, cardiac chair position and one touch zero level (CPR) position operated by separate patient control and nurse control panel.
- Detachable/ Mounting rod polymer moulded head and foot panel
- Provided with 4 sectional polymer moulded safety side rails with separate angle indicators.
- Corner PU buffers absorb impact and reduce damage in handling.
- Provided with heavy duty saline stand with provision for fixing it at all four locations.
- Dual side manual CPR
- Frame size (approx): 2030l x 1060 w x 550-750 h
- Back rest tilt : 0-75 degree
- Knee rest tilt : 0-35 degree
- Trendelenburg : 0-12 degree
- Reverse trendelenburg : 0-12 degree
- The complete metal components are 7 tank pre treated and epoxy powder coated.

Optional: Central Locking System



Feature

s:

- Made up of CRCA tube and CRCA Sheet.
- Available with four section 40 density form mattress
- Castors of 125/150 mm dia. With diagonally locking mechanism.
- Position: Back rest, knee rest, trendelenburg / reverse trendelenburg, high – low, cardiac chair position and one touch zero level (CPR) position operated by separate patient control and nurse control panel.
- Detachable/ Mounting rod polymer moulded head and foot panel
- Saline stand provision on four locations.
- Split type 4 section polymer moulded safety side rails with separate angle indicator
- The complete metal components are pre treated and epoxy powder coated.
- Frame size (Approx): 2030 L * 900 W * 550 – 750 H
- Back rest tilt: 0 - 75 degree
- Knee rest tilt: 0 - 35 degree
- Trendelenburg: 0 - 12 degree
- Reverse Trendelenburg: 0 - 12 degree
- Corner PU buffers absorb impact and reduce damage in handling.



Features:

- Made up of CRCA tube and CRCA Sheet.
- Available with four section 40 density form mattress
- ABS head and foot bows are locked and provide a soothing look.
- Corner PU buffers absorb impact and reduce damage in handling.
- Saline stand provision on four locations.
- Castors of 125mm dia. With diagonal locking are provided for mobility with control.
- The complete metal components are pre treated and epoxy powder coated.
- Back rest, Knee rest, High - Low Positions operated by hand held remote control system.
- Frame size (Approx): 2030 L * 900 W * 550 – 750 H
- Back rest tilt: 0 - 75 degree
- Knee rest tilt: 0 - 35 degree

**• Features:**

- Mattress area divided into 4 sections made up of CRCA Sheet.
- Detachable ABS head and foot bows provide a soothing look.
- 4 pieces individual ABS side rails with angle indicator
- Corner PU buffers absorb impact and reduce damage in handling.
- Saline stand provision on four locations.
- Castors of 125mm dia. With diagonal locking are provided for mobility with control.

- The complete metal components are pre treated and epoxy powder coated.
- Back rest, Knee rest, Trendelenburg/Reverse Trendelenburg, High –Low Positions maneuverable from foot-end by 4 separate crank mechanisms.
- Frame size (Approx): 2030 L * 900 W * 600 – 800 H
- Back rest tilt: 0 - 75 degree
- Knee rest tilt: 0 - 35 degree
- Trendelenburg: 0 - 8 degree
- Reverse Trendelenburg: 0 - 8 degree



Feature

s:

- Steel structure with Mattress area divided into 4 sections provides for 5 different position adjustments.
- Detachable SS framed head and foot bows with laminated panels provide ease of access and a aesthetic look.
- SS/MS Collapsible side railing.
- Corner rubber / PU buffers absorb impact and reduce damage in handling.
- Saline stand provision on four locations.
- Castors of 125mm dia. With diagonal locking are provided for mobility with control.
- The complete metal components are pre treated and epoxy powder coated.
- Positions maneuverable from foot-end by 4 separate crank mechanisms.
- Frame size: 2060 L * 900 W * 600 – 800 H
- Back rest tilt: 0 - 75 degree
- Knee rest tilt: 0 - 35 degree
- Trendlenburg: 0 - 8 degree
- Reverse Trendlenburg: 0 - 8 degree

**Features:**

- Perforated CRC Sheet top in four sections
 - Backrest, Kneerest Positions obtained by electric actuator operator by Hand Remote System
 - ABS Moulded Head & Foot Bows
 - ABS Moulded Safety side railings
 - Telescopic I.V. Rod with 4 locations
 - Bed Mounted on 100mm Deluxe Castors, two with brakes
 - Pre Treated and Epoxy powder Coated.
 - Frame work made of rectangular M.S. tubes
 - Overall approx Size: 2030mm x 900mm x 600 mm
-

FOWLER BED MECHANICAL

MF 004A



Features:

- Mattress area divided into 4 sections made up of CRCA Sheet, mounted on 100 mm castors two with brake.
- Detachable ABS head and foot bows with 4 pieces individual ABS side rails. Head & foot rest adjusted by individual crank.
- Saline stand provision on four locations.
- The complete metal components are pre treated and epoxy powder coated.
- Frame size (Approx): 2030 L * 900 W * 600H
- Back rest tilt: 0 - 75 degree & Knee rest tilt: 0 - 35 degree



FOWLER BED DELUXE

MF 004B

Features:

- Mattress area divided into 4 sections made up of CRCA Bed mounted on 100 mm casters, two with brake.
- Detachable wooden laminated head and foot bows. Head and foot rest adjusted by individual crank
- Saline stand provision on four locations.
- The complete metal components are pre-treated and epoxy powder coated.
- Frame size (Approx): 2060 L * 900 W * 600H
Back rest tilt: 0 - 75 degree & Knee rest tilt: 0

- 35 degree



FOWLER BED STANDARD

MF 004C



Features:

- Mattress area divided into 4 sections.
- Head and foot bows uneven in height with vertical tube supports.
- Saline stand provision on four locations.
- The complete metal components are pre treated and epoxy powder coated.
- 2 separate crank mechanisms for Back rest and Knee rest respectively.
- Back rest tilt: 0 - 75 degree
- Knee rest tilt: 0 - 35 degree

SEMI FOWLER BED

MF 005

Features:

- Mattress area divided into 2 sections.
- Saline stand provision on four locations.
- The complete metal components are pre treated and epoxy powder coated.
- Crank mechanism for Back rest.
- Frame size: 2030 L * 900 W * 600H
- Back rest tilt: 75 degree



SEMI-FOWLER BED DELUXE

MF 005A



Features:

- Mattress area divided into 2 sections.
- Detachable SS framed head and foot bows with laminated panels and back rest adjusted by crank mechanism
- Saline stand provision on four locations.
- The complete metal components are pre treated and epoxy powder coated.
- Frame size: 2030 L * 900 W * 600H
- Back rest tilt: 75 degree



SEMI-FOWLER BED STANDARD

MF 005B

Features:

- Mattress area divided into 2 sections.
- Head and foot bows uneven in height with vertical tube supports.
- Saline stand provision on four locations.
- The complete metal components are pre treated and epoxy powder coated.
- Crank mechanism for Back rest adjustment.
- Frame size: 2030 L * 900 W * 600H
- Back rest tilt: 75 degree



PLAIN WARD BED DELUXE

MF 006



Features:

- Frame made of rectangular CRCA pipes.
- Bed top made of CRCA perforated sheet.
- SS frame head and foot bows with laminated panels.
- Provision for fixing IV/Saline Stand.
- The complete metal components are pre treated and epoxy powder coated.
- Frame size: 2030 L * 900 W * 600H

PLAIN WARD BED STANDARD

MF 006A

Features:

- Frame made of rectangular CRCA pipes.
- Bed top made of CRCA perforated sheet.
- Head and foot bows uneven in height with vertical tube supports provided with mattress guard.
- Provision for fixing IV/Saline Stand.
- The complete metal components are pre treated and epoxy powder coated.

Frame size: 2030 L * 900 W * 600H



ATTENDANT BED

MF 007

Features:

- Frame made of rectangular CRCA pipes.
- Bed top made of CRCA perforated sheet.
- Frame is mounted on PVC stumps.
- The complete metal components are pre treated and epoxy powder coated.

Frame size: 1870 L * 600 W * 450H



ATTENDANT BED CUM CHAIR

MF008



Features:

- Frame is made up of CRCA tubular pipe mounted on PVC stumps.
- Padded with good quality mattress.
- Metal body is pretreated and epoxy powder coated.

Dimension is 1900 * 620 * 580 mm.

PEDIATRIC BED

MF 009

Features:

- Frame made of CRCA rectangular tubes.
- Bed top made of CRCA perforated sheet.
- Head and foot bows made of CRCA pipes with vertical tube supports.
- Drop down/swing down side railing made of CRCA pipes

Frame size: 1370 L * 760 W * 600H Approx.





Features:

- Frame made of rectangular CRCA pipes with SS swing-down railing and Removable radiolucent top.
- 2 section top with ratchet mechanism for back rest.
- Trendlenburg/Reverse Trendlenburg position maneuverable using pneumatic pump.
- Castors of 125mm dia. With diagonal locking provided for mobility with control.
- Oxygen cylinder holder and SS IV stand.
- The complete metal components are pre treated and epoxy powder coated.
- Frame size: 2100 L * 560 W * 650-950mm H

***Also available in complete SS**

***Also available in with High-Low manoeuvre via crank mechanism**

PATIENT TROLLEY- DELUXE

MF 011



Features:

- Removable with Fixed Height
 - Frame is made of CRCA tube.
 - Trolley mounted over 150mm caster for ease of movement.
 - Provision of saline stand.
 - Provided with swing down railing.
 - Provided with cylinder cage
 - Dimension (Approx): 2100 L*560 W*810 H (mm)
- *Also available in complete Stainless steel.**

STRETCHER ON TROLLEY- STANDARD

MF 011A

Features:

- Removable with Fixed Height
- Frame is made of CRCA tube.
- Trolley mounted over 150mm caster for ease of movement.
- Provision of saline stand.
- Dimension (Approx): 2100 L*560 W*810 H (mm)

***Note: Optional mattress**

***Also available in complete Stainless steel.**



**Features:**

- Three sections top made up of SS.
 - Leg section sliding completely under the main section.
 - Base frame made up of CRCA tubes.
 - Trendlenburg's position achieved by pneumatic pump.
 - Back rest adjustable with ratchet mechanism.
 - Lithotomy crutches/Lithotomy rods with straps.
 - Provided with suitable mattress.
 - The complete metal components are pre treated and epoxy powder coated.
 - Dimension [Approx.]:1830 L x 680 W x 760 H mm.
-

**Features:**

- Frame made of CRCA pipes/sheets with adjustable perineal cut.
 - Three sectional top made of Stainless Steel fixed on trolley base mounted on 125 mm castors
 - Leg section sliding completely under the main section.
 - Tendelenburg's position can be achieved by using pneumatic pump.
 - Waste collection tray is positioned below the perineal cut.
 - Back rest adjustable with ratchet mechanism.
 - Height can be adjusted by hydraulic pump.
 - Hand grips, Lithotomy crutches and mattress
 - The complete metal component is pre treated and epoxy powder coated.
 - Dimension [Approx]:1830 L x 680 W x 760 H mm.
-



Features:

- CRCA sheet construction
- Two section Top with Molded Mattress
- Backrest movement on ratchet mechanism
- Five drawers of different size to accommodate different items and one cabinet
- BP apparatus tray
- Leveler at base
- Pre-treated and powder coated

Approx Dimension: 1880 L* 610W *810 H mm

EXAMINATION COUCH- STANDARD

MF 013A



EXAMINATION COUCH- GYNAE

MF 013B

Features:

- Made of CRCA sheets.
- Back rest adjustable with pneumatic pump.
- Padded with foam top in high quality washable Rexene.
- Provided with three drawers and three utility cabinets.
- Provided with foldable BP apparatus tray on the side.
- Provided with inbuilt foot step.
- Provision of saline stand.
- The complete metal component is pre treated and epoxy powder coated.
- Dimension is 1870 (L)* 560 (W) *810 (H) mm



**Features:**

- Mattress area divided into 3 sections made up of CRCA sheet with high density three section mattress
- Detachable polymer moulded head and foot bows provide a soothing look.
- Provided with polymer moulded safety side rails.
- Corner PU buffers absorb impact and reduce damage in handling.
- Provided with heavy duty saline stand with provision for fixing it at all four locations.
- Castors of 125/150 mm dia. With central locking mechanism for mobility with control
- Back rest, trendelenburg / reverse trendelenburg, high – low, position maneuverable by separate actuator operated by patient control.
- Frame size (approx): 2030L X 900 W X 500-750 H (when fully extended)
- Frame size (approx): 2030L X 900 W X 500-750 H (when leg section under the middle section)
- Back rest tilt : 0-75 degree
- Trendelenburg : 0-12 degree & Reverse trendelenburg : 0-12 degree
- The complete metal components are 7 tank pre treated and epoxy powder coated.
- Accessories: one collection tray, one set of lithotomy rod and crutches, one heavy duty iv pole, one set of hand grip

Optional**Dust Protective Base Cover**

EXAMINATION TABLE- DELUXE

MF 015



Features:

- Sheet and pipes are made up of CRCA material.
- Back rest can be adjusted with pneumatic pump.
- Provided with suitable mattress.
- Dimension is 1830 (L)* 610 (W) *810 (H) mm

GYNAE EXAMINATION TABLE DELUXE

MF 015A

Features:

- Two section cushioned top made up of CRCA sheet.
- Frame is made up of CRCA rectangular pipe.
- Back rest adjustable by ratchet mechanism.
- SS waste collection bowl provided below the Perineal cut
- All metal components are pre-treated and epoxy powder coated.
- Dimension Approx: 1830 L* 610W *810 H mm.



OVERBED TABLE WOODEN TOP

MF 016



Features:

- Rectangular framework mounted on 5 cm castors.
- Wooden top adjustable by Gear.
- All metal components are pre treated and epoxy powder coated.
- Dimension: 760mm L * 400mm
W*760mm-1050mm

Optional: Available with Pneumatic Pumps (DELUXE)

Available with Knob (STANDARD)

OVERBED TABLE ABS TOP

MF 016A

Features:

- Overbed table with ABS top.
- Size: 770 × 380 × 750 / 1010 mm.
- Framework is made of Epoxy steel.
- Adjustable height from 750 - 1010 mm by pneumatic pump
- 4 Swivel castors all with brakes.

**** Optional: Gear Handle**



BEDSIDE LOCKER DELUXE

MF 017



Features:

- Made up of CRCA pipe/sheets with SS top.
- Provided with one storage box and one drawer.
- CRCA tubular leg, two is mounted on PVC stump and two on 5cms castors.
- Pre-treated and epoxy powder coated.
- Approx Dimension: 405mm L* 405mm W* 810mm H.

MEDICINE SIDE CUPBOARD- STANDARD

MF 017A



Features:

- Made up of CRCA pipe/sheets with SS top.
- Provided with one storage box.
- CRCA tubular leg Mounted on PVC stump.
- Pre-treated and epoxy powder coated.
- Approx. Dimension: 405mm L* 405mm W* 810mm H.

MEDICINE SIDE CUPBOARD- ULTRA PLUS

MF 017B



Features:

- Made up of CRCA pipe/sheets with SS top.
- Provided with two shelves and two drawers.
- Mounted on 5 cm casters.
- Pre-treated and epoxy powder coated.
- Approx. Dimension: 405mm
- L* 405 mm W* 810 mm H.

MEDICINE SIDE CUPBOARD

MF 017C

Features:

- Made up of CRCA pipe/sheets.
- Top / Drawer / cabinet is made of laminated / polymer moulded material.
- One drawer & one cabinet with handle facility.
- Mounted on four 50 mm PVC castors.
- Complete metal component is pre-treated & epoxy powder coated.
- Dimension: 405mm L* 405 mm W* 4810 mm H.



Baby Cradle - Swinging Type

MF 018



Features:

- Size: L 910L x W 380 x H 990mm.
- Stand made of ERW tubular pipe with hanging hooks mounted on rubber shoes.
- CRIB made of steel bars/strips.
- Pre-treated and epoxy powder coated.

CRASH CART FULLY SS / MS CRASH CART FULLY SS / MS

MF 019



Features:

- Framework is made up of SS Tubes.
 - SS storage cabinet with three drawers.
 - Provision of M.S oxygen cylinder holder.
 - Six colour removable bins.
 - Provided with three S.S shelves.
 - Provided with S.S I.V rod and corner buffer and cardiac massage board.
 - The cart is mounted on 125 mm dia castors 2 with brakes.
-

CRASH CART MS/PLASTIC CABINET

MF 019A



- Framework is made up of CRCA tubes.
- Provided with two polystyrene boxes each with three lockable drawers.
- Provision of M.S oxygen cylinder holder.
- Six colour removable bins.
- Provided with three S.S Shelves.
- Provided with S.S I.V rod and corner buffer and cardiac massage board.
- The cart is mounted on 125 mm dia castors 2 with brakes.
- Approx. Dimension: 960 mm L * 500 mm W * 1545 mm H.

**** Also available in S.S framework and PC Cabinet.

MONITOR TROLLEY PC

MF 020



Features:

- Made up of CRCA sheet with four shelves and one drawer.
- Mounted on 100 mm dia castors two with break.
- Trolley has socket to receive 5/15 Amp – 230V A/C 3 – pin plugs. Two inputs each with safety fuse and switch.
- Trolley has two handle for ease of relocation.
- All metal components are pre treated and epoxy powder coated.
- Approx Dimension: 680 mm L * 610W mm * 1500 mm H

*** Also available in all Stainless Steel.

MEDICINE TROLLEY

MF 021

Features:

- Framework is made up of CRCA tubes mounted on 10 cm castors.
- Two SS shelves with three sides railing on top self and 4 side railing in lower self.
- Two drawers under each self.
- All metal components are pre-treated and epoxy powder coated.
- Dimension: 760 mm L * 460 mm W* 810 mm H.

*Also available in all Stainless Steel.



DRESSING TROLLEY

MF 022



Features:

- Made up of SS tubes mounted on 100 mm castors.
- Two SS shelves top with 3 side railing and bottom with 4 side railing.
- Bowl and bucket holding cage on one side with SS bowl and bucket.
- Dimension: 1010 mm L x 460 mm W x 810 mm H.

INSTRUMENT TROLLEY

MF 023

Features:

- Made up of SS pipes and sheets mounted on 100 mm castors.
- Two SS shelves, top self with 3 side railing and lower self with 4 side railing.
- Dimension:

(MF023A) 600 mm L x 460 mm W x 810 mm H

(MF023B) 760 mm L x 460 mm W x 810 mm H

(MF023C) 910 mm L x 460 mm W x 810 mm H



MAYO'S TROLLEY WITH GEAR

MF 024



Features:

- Made up of CRCA pipes mounted on 5 cms castors.
- Height adjustable by gear provided on side.
- Dimension of SS tray: 600 mm L x 460 mm W.

MAYO'S TROLLEY DOUBLE STAND WITH GEAR

MF 024A

Features:

- Framework is made up of SS Tubes mounted on 100 mm castors.
- Height adjustable by gear provided on either side.
- Dimension of SS tray: 910 mm L x 500 mm W.



I.V. STAND

MF 025



- MF 25A- SS pipe with four leg SS base.
- MF 25B- SS pipe with plastic base with 4 hooks.
- MF 25C- MS pipe with 5 legs Plastic base.

CYLINDER TROLLEY

MF 026



Features:

- Made up of CRCA pipes mounted.
- Pre-treated and epoxy powder coated.
-
- A. For D type cylinder mounted on 100mm castors.
- B. For B type cylinder mounted on 150mm castors.

BASIN STAND DOUBLE



MF 027A

Features:

- Made up of CRCA pipes mounted on 10 cms castors.
- Provided with 2 SS basin.
- Dimension: 350mm dia * 850mm H.

*****MF 027B Also available in SS frame.**

BASIN STAND SINGLE

Features:

- Made up of CRCA pipes.
- Five leg PVC base mounted on 5cms dia.
- Provided with one SS basin.
- Dimension: 350mm dia * 850mm H.

*****MF 027D Also available in SS.**

MF 027C



KICK BUCKET



MF 027E

Features:

- SS Bucket with stainless steel frame with castors.
- Dimension: 35cm D x 35cm H

FOOT STEP DOUBLE SS

MF 028A



Features:

- Frame made up of SS tube fitted on PVC stumps.
- Step made up of SS Sheet.
- Dimension: 18" Lx18"W x18"H.

First step height 230mm

Second step height 450mm

FOOT STEP SINGLE SS

MF 028B

Features:

- Made from SS Shuts mounted on PVC stumps.
- Also available in CRCA Shuts Epoxy-Powder Coated.
- Dimension: 18" Lx 9"W x 9"H.



FOOT STEP DOUBLE MS

MF 028C



Features:

- Made from CRCA Sheets with rubber sheet on top mounted on PVC Stumps.
- Pre treated and Epoxy Powder coated.
- Dimension: 18" Lx 18"W x18"H.
- First step height 230mm
- Second step height 450mm

***MF 028D Also available in complete SS**

SURGEON'S STOOL DELUXE

MF 029A



Features:

- Cushioned top with 5 leg based mounted on 5 cms castors.
- Adjustable height with pneumatic pump.
- Provided with back support & tubular ring for foot support.

SURGEON'S STOOL GENERAL

MF 029



MF 029B



MF 029C

Features:

- Cushioned top with 5 leg based mounted on 5 cms castors.
- Adjustable height with pneumatic pump.
- Tubular ring for foot support.

*****Also available in MS frame.**

INSTRUMENT CABINET (MS)

MF 031



Features:

- Made of CRCA sheets provided with single door and glass shelves.
- Pre-treated and epoxy powder coated.
- A - 63"Hx24"Wx18"D WITH 5 GLASS SHELVES(6 CHAMBERS)
- B - 58"Hx24"Wx14"D WITH 4 GLASS SHELVES(5 CHAMBERS)

BEDSIDE SCREEN (4 FOLDS WITH CURTAINS)

MF 032



Features:

- Made of CRCA pipes on 5 cms castors provided with-out curtains/springs. Powder coated.
- Dimension: 66" H x 96" L.

***Option: MF 032A Also available in 3 fold**

SOILED LINEN TROLLEY (CANVAS BAG)

MF 033



Features:

- SS framework mounted on 10 cms castors.
- Provided with superior quality canvas bag.

**** OPTION : MF 033A Also available in CRCA frame work**

SOILED LINEN TROLLEY (PLASTIC BUCKET)

MF 033A

Features:

- SS framework mounted on 10 cms castors.
- Provided with superior quality plastic bucket.
- Push handle for ease of handling.

**** MF 033B Also available in CRCA framework**



LINEN TROLLEY WITH THREE SHELVES AND BAG

MF 034A



Features:

- Made of CRCA pipes/sheet
- Provided with three shelves and one canvas bag with diameter of 460 mm
- Mounted on 100 mm castors.
- The complete metal components are pre treated and epoxy powder coated.
- Approx. Dimensions: L 960 mm x W 500 mm x H 1060

LINEN TROLLEY WITH CANVAS BAGS

MF 034B

Features:

- Made of SS/MS Round Pipes.
- Provided with canvas bags.
- Mounted on 100 mm castors.
- Approx. Dimensions: L 960 mm x W 500 mm x H 910 mm

**** Optional: 2 Canvas Bag**



BIO MEDICAL WASTE BINS

MF 035



Features:

- 22 LTRS
- Colour: Red, Blue, Black, Green, Yellow, White
- With S.S. Trolley
- Material: Virgin Polymers
- Applications: Hospitals, Nursing homes, Clinics, Dental Hospital, ENT Hospitals.

STRETCHERS

MF 036



- MF 036B Folding Stretcher 2 Fold

- MF 036A Folding Stretcher 4 fold(Canvas)



SCOOP STRETCHER

MF 037



WHEELCHAIR NON-FOLDING

MF 038



MF 038A



MF038B

Features:

- Non folding 17" seat and fixed cushioned back rest, fixed arm-rest, 6" front castors, 26" rear wheel with brakes.
 - A) SS Dimension: 38"Lx25"Wx36" H
 - B) MS Pre-treated epoxy powder coated.
 - Dimension: 38"Lx25"Wx36" H
-

WHEELCHAIR FOLDING

MF 038C



Features:

- Foldable wheelchair 17" square seat, fixed arm-rest, 26" rear wheel, 6" front castors.
- Pre-treated powder coated.
- Dimension: 42" L x 26" W x 36" H.

THREE SEATER WAITING CHAIR

MF 039



Features:

- Perforated Sheet
- Chrome Plated Finish
- Size (Approx)- 69"Lx 25"W x 30" H

***Also Available in SS**

SANARE MEDICAL



Website :- www.sanaremedical.com

Email :- info@sanaremedical.com

Tel :- +91 9810275300, +91 87500 57290

Corp.Office :- E-354 Greater Kailash Part-1 ,New Delhi-110048

Head Office :- DP 61, Phase 1 Udyog Vihar Gurugram Haryana-122001