




**SANARE
MEDICAL**



PRODUCT CATALOG



+91 011-41024178, 011-41437486



www.sanaremedical.com



Mobility Aids

SINGLE STICK

Product Code : MESTK01


ELBOW CRUTCHES

Product Code : MEELBC1


TRIPOD STICK


Product Code : METSK01


QUAD STICK

Product Code : MEQSK01


UNDERARM CRUTCHES

Product Code : MEUCC01


Features

- Height adjustable
- Chrome / Aluminium Frame

- Anti-slipping Bushes
- Superior Quality

WALKER REGULAR (POWDER COATED)

Product Code : MQFWR01



Features

- Walker with reciprocal feature
- Compact storage & lateral access
- One-button for easy folding
- 6 pin height adjustment
- Tool free adjustments

Walker Imported (Chrome finish)

Product Code : MEFWR01



- Chrome frame
- Tool free adjustments
- One-button for easy folding

Walker with Wheels

Product Code : MEFWW01



Features

- Height Adjustments
- Compact Storage
- Lateral Access

MULTI FUNCTION WALKER

Product Code : MEFWR03


**6 in 1
product**


Lifting Aid



Walker



Shower Chair



Commode Chair



Toilet Safety Frame



Toilet Seat Raiser

ROLLATOR CUM WHEELCHAIR

Product Code : MERTB02


Features

- Powder Coated Frame
- Bag or Basket option
- 8" castors
- Loop style handle brake (push down to brake)
- Foldable and can fit into the car easily
- Removable padded backrest

ROLLATOR WITH FOOTREST

Product Code : MERTB02


3 in 1 Wheelchair


Chair

Transport Wheelchair

Walker

for short distances

Height Adjustable Commode Chair Powder Coated

Product Code : MEPHAC1



Height Adjustable Commode Chair with Wheels (Powder Coated)

Product Code : MEPCWC1



Height Adjustable Commode Chair (Chrome)

Product Code : MEFCH01



Height Adjustable Commode Chair with Wheels (Chrome)

Product Code : MEFCW01



Features

- Easily opens and folds
- Folds flat for convenient storage
- Durable plastic snap on seat and lid
- Plastic arms for added comfort
- Comes complete with commode bucket, cover and is rust resistant
- 2 rear locking(optional) swivel casters

**Commode Stool with lock
Powder Coated**

Product Code : MEFCL01



Pot extra

**Commode Chair with
armrest Powder Coated**

Product Code : MEFCA01


OVERBED FOOD TABLE

Product Code : MEFTB01


Features

- Easy to adjust height
- Long lasting material
- Deluxe wheels 4 side lock system
- Wooden frame for long life
- Non Rusting Powder Coated
- Light-weight
- Superb finish Application

PREMIUM BASIC WHEELCHAIR CHROME

Product Code : MEPBBC1





Features

- Premium Steel Frame. Chrome Finish
- Premium Washable Upholstery
- Fixed Armrest
- Height Adjustable Footrest
- Solid Composite Castors
- Solid **Rubber** Rear Wheels
- Dual Brakes
- Fixed Seatbelt & calf strap
- Chrome finish
- Available in 2 colors

PREMIUM BASIC WHEELCHAIR POWDER COATED

Product Code : MEPRWC2



Features

- Black Powder Coated Frame
- Premium Nylon Upholstery
- Fixed PVC padded armrest
- Height Adjustable Footrest
- 8" front Solid Composite Castors
- 21" rear spoke and aluminium Solid Wheels
- Dual Brakes
- Fixed Seatbelt
- Durable Plastic footrest

Made in India



ASTM Standards



1 Year Warranty*
Terms & Conditions apply

WHEELCHAIR WITH U CUT COMMODE

Product Code : MECWC01



Made in India



ASTM Standards


 Year Warranty*
Terms & Conditions apply

Features

- Steel Chrome / Powder coated Frame
- PVC Upholstery
- Fixed PVC padded Armrest
- Durable Plastic footrest
- Footrest with calf strap
- 8" front Solid Composite Castors
- 24" rear spoke and aluminium Solid Wheels
- Dual Brakes

WHEELCHAIR WITH SEAT LIFT COMMODE

Product Code : MECWC02



Made in India



ASTM Standards


 Year Warranty*
Terms & Conditions apply

Features

- Steel Chrome / Powder coated Frame
- PVC Upholstery
- Fixed PVC padded Armrest
- Durable Plastic footrest
- Footrest with calf strap
- 8" front Solid Composite Castors
- 24" rear spoke and aluminium Solid Wheels
- Dual Brakes

RECLINING WHEELCHAIR WITH COMMODE (U CUT SEAT)

Product Code : MERWC01



U cut seat

Features

- Premium Steel Frame
- Detachable PVC padded Armrest
- Deluxe U cut toilet seat
- Detachable and elevating Footrest
- 8" front wheels with solid tyre
- 24" rear spoke and aluminium wheels
- Aluminium footplate & chrome handrim
- Premium hydraulic/
Slotted reclining system

RECLINING WHEELCHAIR WITH COMMODE (SEAT LIFT)

Product Code : MERWC02



Seat lift



Seat lift
commode

Features

- Premium Steel Frame
- Detachable PVC padded Armrest
- Deluxe U cut toilet seat
- Detachable and elevating Footrest
- 8" front wheels with solid tyre
- 24" rear spoke and aluminium wheels
- Aluminium footplate & chrome handrim
- Premium hydraulic/
Slotted reclining system

TRANSIT WHEELCHAIR

Product Code : METWC05



Features

- Premium Aluminum Frame
- Flipback back rest
- Fixed foam padded detachable Armrest
- Flipup Footrest
- 6" front PVC castor wheels
- 12"/14" solid mag rear wheels
- Aluminium footplate
- Premium nylon upholstery -black/checks
- Weight bearing capacity 100 kilos*
- Fixed seat belt

BARIATRIC WHEEL CHAIR-CHROME

Product Code : MOBWC01



Cushion extra

Features

- Heavy duty chrome frame with Dual chassis
- Seat width 22"
- Detachable and elevating footrest
- Heavy duty PU front castor wheel
- 24" rear spoke and aluminum wheel with black handrim
- Optional Premium Cushion for maximum comfort
- Weight bearing capacity 150 kilos*
- Wheelchair designed for big and tall

Deluxe Wheelchair with PU Mag Wheels

Product Code : MEBWM01



Features

- Premium Powder Coated / Chrome Frame
- Premium Washable Upholstery
- Fixed PVC padded Armrest
- 8" front castor solid tyre
- 24" rear Mag solid figured tyre
- Solid Rubber Rear Wheels
- Dual Brakes
- Fixed Seatbelt

PEDIATRIC WHEELCHAIR

Product Code : MEPWC01



Features

- Solid powder coated steel frame
- Seat width 14"
- Specially designed for children
- Fixed black PVC padded armrest
- Fixed aluminum footplate
- 6" PVC front castor
- 20" rear solid spoke wheel
- Maximum weight capacity 80kilos

BASIC ELECTRIC WHEEL CHAIR

Product Code : MEEWC01



Technical Specifications

Driving range : 26-31 km (depending on user weight and road condition)	Loading capacity: 100kg
Max speed : 6km/h	Front wheel : 8 inch
N.W. : 36kgs	Rear wheel : 12 inch
G.W. : 42kgs	Motor : Spur Gear Motor
Weight without battery: 32kgs	Battery : 12V 20Ah
Overall size : 1170*620*940mm	Controller : Intelligent control
Folded size : 745*370*730mm	Climbing ability : 8°
Seat Width : 450(L)*430(W)mm	Brake : Hand brake
	Warranty : 1 year

RECLINING ELECTRIC WHEEL CHAIR

Product Code : MEWCR01



Technical Specifications

Driving range : 26-31km (depending on user weight and road condition)	Front wheel : 8 inch
Max speed : 6km/h	Rear wheel : 22 inch
N.W. : 49kgs	Motor : Spur Gear Motor
G.W. : 55kgs	Battery : 12V 20Ah
Weight without battery: 45kgs	Controller : Intelligent control
Overall size : 1090*685*1275mm	Climbing ability : 8°
Folded size : 805*420*730mm	Brake : Hand brake
Seat Width : 450(L)*430(W)mm	Frame color options : Silver
Loading capacity : 100kg	Warranty : 1 year



Head Office:

**DP-61 Phase 1 Udyog Vihar
Gurugram Haryana 122001 INDIA**



+91 98102 75300 / 87500 57290



info@sanaremedical.com



www.sanaremedical.com